

Product Datasheet

Bovine MCP2 ELISA Kit (orb1197126)

Description	Bovine Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q09141
Sensitivity	6.4pg/mL
Entrez	281044
Expiration Date	Please enquire.