

Product Datasheet

Porcine MCP1 ELISA Kit (orb1197127)

Description	Porcine Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit
Reactivity	Porcine
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42831
Sensitivity	0.124ng/mL
Entrez	397422
Expiration Date	Please enquire.