

Product Datasheet

Porcine LTF ELISA Kit (orb1197128)

| | |
|----------------------------|---|
| Description | Porcine Lactoferrin (LTF) ELISA Kit |
| Reactivity | Porcine |
| Range | 31.2-2000pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma or other biological fluids. |
| Uniprot ID | P14632 |
| Sensitivity | 13.2pg/mL |
| Entrez | 397649 |
| Expiration Date | Please enquire. |