

Product Datasheet

Bovine LEP ELISA Kit (orb1197135)

Description	Bovine Leptin (LEP) ELISA Kit
Reactivity	Bovine
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P50595
Sensitivity	0.063ng/mL
Entrez	280836
Expiration Date	Please enquire.