

## Product Datasheet

### Porcine HIF1a ELISA Kit (orb1197149)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Porcine Hypoxia Inducible Factor 1 Alpha (HIF1a) ELISA Kit    |
| <b>Reactivity</b>          | Porcine   |
| <b>Range</b>               | 15.6-1000pg/mL  |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | sandwich  |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <a href="#">A0T1H6</a>  |
| <b>Sensitivity</b>         | 5.8pg/mL  |
| <b>Entrez</b>              | <a href="#">396696</a>  |
| <b>Expiration Date</b>     | Please enquire.   |