

## Product Datasheet

### Bovine CD40L ELISA Kit (orb1197194)

<b>Description</b>	Bovine Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	0.625-40ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P51749</a>
<b>Sensitivity</b>	0.18ng/mL
<b>Entrez</b>	<a href="#">282387</a>
<b>Expiration Date</b>	Please enquire.