

Product Datasheet

Bovine bEP ELISA Kit (orb1197207)

Description	Bovine Beta-Endorphin (bEP) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01190
Sensitivity	4.89pg/mL
Entrez	281416
Expiration Date	Please enquire.