

Product Datasheet

Bovine AP ELISA Kit (orb1197213)

Description	Bovine Aprotinin (AP) ELISA Kit
Reactivity	Bovine
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P00974
Sensitivity	0.43ng/mL
Entrez	404172
Expiration Date	Please enquire.