

Product Datasheet

Bovine ANXA5 ELISA Kit (orb1197214)

Description	Bovine Annexin A5 (ANXA5) ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P81287
Sensitivity	0.16ng/mL
Entrez	281626
Expiration Date	Please enquire.