

Product Datasheet

Bovine ITaC ELISA Kit (orb1197236)

Description	Bovine Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
Reactivity	Bovine
Range	46.8-3000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	A9QWQ1
Sensitivity	15.46pg/mL
Entrez	516104
Expiration Date	Please enquire.