

Product Datasheet

Bovine IL22 ELISA Kit (orb1197237)

Description	Bovine Interleukin 22 (IL22) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	A5JS75
Sensitivity	5.4pg/mL
Entrez	507778
Expiration Date	Please enquire.