

## Product Datasheet

### Rat NT ELISA Kit (orb1197294)

<b>Description</b>	Rat Neurotensin (NT) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	4.68-300pmol/L
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>P20068</b>
<b>Sensitivity</b>	1.3pmol/L
<b>Entrez</b>	<b>299757</b>
<b>Expiration Date</b>	Please enquire.