

Product Datasheet

Rat APOC2 ELISA Kit (orb1197422)

Description	Rat Apolipoprotein C2 (APOC2) ELISA Kit
Reactivity	Rat
Range	39-2500pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	G3V8D4
Sensitivity	14pg/mL
Entrez	292697
Expiration Date	Please enquire.