

## Product Datasheet

### Mouse UCHL1 ELISA Kit (orb1197482)

<b>Description</b>	Mouse Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	78.1-5000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>Q9R0P9</b>
<b>Sensitivity</b>	25pg/mL
<b>Entrez</b>	<b>22223</b>
<b>Expiration Date</b>	Please enquire.