

Product Datasheet

Mouse TRPM4 ELISA Kit (orb1197495)

Description	Mouse Transient Receptor Potential Cation Channel Subfamily M, Member 4 (TRPM4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q7TN37
Sensitivity	0.048ng/mL
Entrez	68667
Expiration Date	Please enquire.