

Product Datasheet

Mouse SUOX ELISA Kit (orb1197529)

Description	Mouse Sulfite Oxidase (SUOX) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8R086
Sensitivity	0.054ng/mL
Entrez	211389
Expiration Date	Please enquire.