

Product Datasheet

Mouse STK11 ELISA Kit (orb1197533)

Description	Mouse Serine/Threonine Kinase 11 (STK11) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9WTK7
Sensitivity	0.051ng/mL
Entrez	20869
Expiration Date	Please enquire.