

Product Datasheet

Canine SPA ELISA Kit (orb1197544)

Description	Canine Surfactant Associated Protein A (SPA) ELISA Kit
Reactivity	Canine
Range	1.87-120ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates or other biological fluids.
Uniprot ID	Q6RXL1
Sensitivity	0.62ng/mL
Entrez	407213
Expiration Date	Please enquire.