

Product Datasheet

Mouse SFRP2 ELISA Kit (orb1197572)

Description	Mouse Secreted Frizzled Related Protein 2 (SFRP2) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P97299
Sensitivity	0.237ng/mL
Entrez	20319
Expiration Date	Please enquire.