

Product Datasheet

Mouse SFRP1 ELISA Kit (orb1197573)

Description	Mouse Secreted Frizzled Related Protein 1 (SFRP1) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q8C4U3
Sensitivity	13.3pg/mL
Entrez	20377
Expiration Date	Please enquire.