

Product Datasheet

Mouse PDIA4 ELISA Kit (orb1197641)

| | |
|----------------------------|---|
| Description | Mouse Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit |
| Reactivity | Mouse |
| Range | 0.312-20ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates or other biological fluids. |
| Uniprot ID | P08003 |
| Sensitivity | 0.112ng/mL |
| Entrez | 12304 |
| Expiration Date | Please enquire. |