

Product Datasheet

Canine NOS2 ELISA Kit (orb1197672)

Description	Canine Nitric Oxide Synthase 2, Inducible (NOS2) ELISA Kit
Reactivity	Canine
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O62699
Sensitivity	0.126ng/mL
Entrez	403822
Expiration Date	Please enquire.