

Product Datasheet

Canine MCP1 ELISA Kit (orb1197715)

Description	Canine Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit
Reactivity	Canine
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P52203
Sensitivity	0.058ng/mL
Entrez	403981
Expiration Date	Please enquire.