

Product Datasheet

Canine LIF ELISA Kit (orb1197735)

Description	Canine Leukemia Inhibitory Factor (LIF) ELISA Kit
Reactivity	Canine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	J9NRL6
Sensitivity	14.7pg/mL
Entrez	403449
Expiration Date	Please enquire.