

## Product Datasheet

### Canine HMG1 ELISA Kit (orb1197790)

<b>Description</b>	Canine High Mobility Group Protein 1 (HMG1) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	3.12-200ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q6YKA4</a>
<b>Sensitivity</b>	1.24ng/mL
<b>Entrez</b>	<a href="#">403170</a>
<b>Expiration Date</b>	Please enquire.