

Product Datasheet

Mouse HADHb ELISA Kit (orb1197799)

Description	Mouse Hydroxyacyl Coenzyme A Dehydrogenase Beta (HADHb) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q99JY0
Sensitivity	0.115ng/mL
Entrez	231086
Expiration Date	Please enquire.