

Product Datasheet

Guinea ECF ELISA Kit (orb1197875)

Description	Guinea pig Eosinophil Chemotactic Factor (ECF) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	G8XRI2
Sensitivity	6.2pg/mL
Entrez	100038010
Expiration Date	Please enquire.