

Product Datasheet

Mouse CYP27B1 ELISA Kit (orb1197909)

Description	Mouse Cytochrome P450 27B1 (CYP27B1) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O35084
Sensitivity	0.214ng/mL
Entrez	13115
Expiration Date	Please enquire.