

## Product Datasheet

### Canine CAV1 ELISA Kit (orb1197972)

<b>Description</b>	Canine Caveolin 1 (CAV1) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	0.156-10ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P33724</a>
<b>Sensitivity</b>	0.064ng/mL
<b>Entrez</b>	<a href="#">403980</a>
<b>Expiration Date</b>	Please enquire.