

## Product Datasheet

### Mouse ART1 ELISA Kit (orb1198012)

<b>Description</b>	Mouse ADP Ribosyltransferase 1 (ART1) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.781-50ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q60935</a>
<b>Sensitivity</b>	0.32ng/mL
<b>Entrez</b>	<a href="#">11870</a>
<b>Expiration Date</b>	Please enquire.