

## Product Datasheet

### Mouse APLN ELISA Kit (orb1198022)

<b>Description</b>	Mouse Apelin (APLN) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	125-8000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<b>Q9R0R4</b>
<b>Sensitivity</b>	37.5pg/mL
<b>Entrez</b>	<b>30878</b>
<b>Expiration Date</b>	Please enquire.