

Product Datasheet

Monkey IFN β ELISA Kit (orb1198066)

Description	Monkey Interferon Beta (IFN β) ELISA Kit
Reactivity	Monkey
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	B6CK14
Sensitivity	7.1pg/mL
Entrez	711414
Expiration Date	Please enquire.