

## Product Datasheet

### Human WISP2 ELISA Kit (orb1198127)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human WNT1 Inducible Signaling Pathway Protein 2 (WISP2) ELISA Kit |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 0.312-20ng/mL  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | sandwich   |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.                     |
| <b>Uniprot ID</b>          | <b>076076</b>  |
| <b>Sensitivity</b>         | 0.133ng/mL   |
| <b>Entrez</b>              | <b>8839</b>  |
| <b>Expiration Date</b>     | Please enquire.  |