

Product Datasheet

Human VCC1 ELISA Kit (orb1198142)

Description	Human VEGF Co Regulated Chemokine 1 (VCC1) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6UXB2
Sensitivity	4.7pg/mL
Entrez	284340
Expiration Date	Please enquire.