

Product Datasheet

Human TRAF3 ELISA Kit (orb1198197)

Description	Human TNF Receptor Associated Factor 3 (TRAF3) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13114
Sensitivity	29pg/mL
Entrez	7187
Expiration Date	Please enquire.