

Product Datasheet

Human SOSTDC1 ELISA Kit (orb1198316)

| | |
|----------------------------|--|
| Description | Human Sclerostin Domain Containing Protein 1 (SOSTDC1) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates or other biological fluids. |
| Uniprot ID | Q6X4U4 |
| Sensitivity | 0.059ng/mL |
| Entrez | 25928 |
| Expiration Date | Please enquire. |