

Product Datasheet

Human SIRT6 ELISA Kit (orb1198352)

Description	Human Sirtuin 6 (SIRT6) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8N6T7
Sensitivity	30.6pg/mL
Entrez	51548
Expiration Date	Please enquire.