

## Product Datasheet

### Human PLP2 ELISA Kit (orb1198539)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Proteolipid Protein 2, Colonic Epithelium Enriched (PLP2) ELISA Kit |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 0.156-10ng/mL   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | sandwich  |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.                            |
| <b>Uniprot ID</b>          | <a href="#">Q04941</a>  |
| <b>Sensitivity</b>         | 0.065ng/mL  |
| <b>Entrez</b>              | <a href="#">5355</a>  |
| <b>Expiration Date</b>     | Please enquire.   |