

## Product Datasheet

### Human PFN1 ELISA Kit (orb1198570)

<b>Description</b>	Human Profilin 1 (PFN1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.1-5000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<b>P07737</b>
<b>Sensitivity</b>	31pg/mL
<b>Entrez</b>	<b>5216</b>
<b>Expiration Date</b>	Please enquire.