

## Product Datasheet

### Human PDIA4 ELISA Kit (orb1198579)

<b>Description</b>	Human Protein Disulfide Isomerase A4 (PDIA4) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>P13667</b>
<b>Sensitivity</b>	0.66ng/mL
<b>Entrez</b>	<b>9601</b>
<b>Expiration Date</b>	Please enquire.