

Product Datasheet

Human PAOX ELISA Kit (orb1198602)

| | |
|----------------------------|---|
| Description | Human Polyamine Oxidase (PAOX) ELISA Kit |
| Reactivity | Human |
| Range | 78.1-5000pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates or other biological fluids. |
| Uniprot ID | Q6QHF9 |
| Sensitivity | 31.5pg/mL |
| Entrez | 196743 |
| Expiration Date | Please enquire. |