

## Product Datasheet

### Human NPSR ELISA Kit (orb1198669)

|                            |                                                |
|----------------------------|------------------------------------------------|
| <b>Description</b>         | Human Neuropeptide S Receptor (NPSR) ELISA Kit |
| <b>Reactivity</b>          | Human                                          |
| <b>Range</b>               | 0.312-20ng/mL                                  |
| <b>Tested Applications</b> | ELISA                                          |
| <b>Note</b>                | For research use only                          |
| <b>Application notes</b>   | sandwich                                       |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>Q6W5P4</b>                                  |
| <b>Sensitivity</b>         | 0.108ng/mL                                     |
| <b>Entrez</b>              | <b>387129</b>                                  |
| <b>Expiration Date</b>     | Please enquire.                                |