

## Product Datasheet

### Human NDUFS1 ELISA Kit (orb1198705)

<b>Description</b>	Human NADH Dehydrogenase Ubiquinone Fe-S Protein 1 (NDUFS1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P28331</a>
<b>Sensitivity</b>	0.105ng/mL
<b>Entrez</b>	<a href="#">4719</a>
<b>Expiration Date</b>	Please enquire.