

Product Datasheet

Human LRP5 ELISA Kit (orb1198813)

Description	Human Low Density Lipoprotein Receptor Related Protein 5 (LRP5) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O75197
Sensitivity	0.118ng/mL
Entrez	4041
Expiration Date	Please enquire.