

## Product Datasheet

### Human HSPG2 ELISA Kit (orb1198930)

<b>Description</b>	Human Heparan Sulfate Proteoglycan 2 (HSPG2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	125-8000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>P98160</b>
<b>Sensitivity</b>	49pg/mL
<b>Entrez</b>	<b>3339</b>
<b>Expiration Date</b>	Please enquire.