

Product Datasheet

Human HSPBP1 ELISA Kit (orb1198932)

Description	Human Heat Shock 70kDa Binding Protein 1 (HSPBP1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9NZL4
Sensitivity	0.105ng/mL
Entrez	23640
Expiration Date	Please enquire.