

Product Datasheet

Human HSP90aA1 ELISA Kit (orb1198939)

Description	Human Heat Shock Protein 90kDa Alpha A1 (HSP90aA1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07900
Sensitivity	0.55ng/mL
Entrez	3320
Expiration Date	Please enquire.