

Product Datasheet

Human HMGA2 ELISA Kit (orb1198960)

| | |
|----------------------------|---|
| Description | Human High Mobility Group AT Hook Protein 2 (HMGA2) ELISA Kit |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates, cell lysates or other biological fluids. |
| Uniprot ID | P52926 |
| Sensitivity | 0.114ng/mL |
| Entrez | 8091 |
| Expiration Date | Please enquire. |