

## Product Datasheet

### Human HBb ELISA Kit (orb1198979)

<b>Description</b>	Human Hemoglobin Beta (HBb) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100µg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or erythrocyte lysates.
<b>Uniprot ID</b>	<a href="#">P68871</a>
<b>Sensitivity</b>	0.45µg/mL
<b>Entrez</b>	<a href="#">3043</a>
<b>Expiration Date</b>	Please enquire.