

Product Datasheet

Human HADHb ELISA Kit (orb1198984)

| | |
|----------------------------|---|
| Description | Human Hydroxyacyl Coenzyme A Dehydrogenase Beta (HADHb) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates or other biological fluids. |
| Uniprot ID | P55084 |
| Sensitivity | 0.053ng/mL |
| Entrez | 3032 |
| Expiration Date | Please enquire. |